




Standard Reference Materials

Increase user-efficiency without compromising on quality

The Standard Reference provide both calibration and verification options for laboratory test equipment, all from the same sample, with their values certified to the highest level of accreditation integrity available globally. The multiple test parameters provide a flexible and cost-effective solution for any laboratory.

Viscosity Standard Oil



Product Description

- Viscosity Standards Oil are for the calibration and verification of viscometers and other viscosity measuring equipment where operating temperature is controlled precisely.
- Fully traceable to IS unit

Packaging type	Plastic Bottle
Packaging size	500mL / bottle
Purpose	calibration/verification
Form	Liquid
Price	Contact Ascendent

Required more details and quotation?

**Please contact sales representative
Ms. Joyce @ +6012-4173813**